

ANTI-OVARY ANTIBODY TEST SYSTEM

For In Vitro Diagnostic Use

AD OMO48 48 Tests
AD OMO96 96 Tests

Introduction:

Ovarian antibodies are found in patients with premature ovarian failure. These patients have an increased frequency of autoimmune diseases; organ and non-organ specific, thyroid being the most common. Women undergoing in-vitro fertilization (IVF) can produce IgG, IgM & IgA antibodies to the theca interna and atretic follicles. Antibodies to the corpora lutea are found in patients with primary sterility and endometriosis.

Principles:

The primary reaction involves circulating antibodies present in the patient's serum which attach to their homologous antigens. This occurs during the incubation period while the serum covers the antigen surface. A secondary reaction then follows in rinsing period which removes all unbound human antibody. The antigen surface is then thoroughly rinsed free of unbound conjugate and viewed under an appropriate fluorescence microscope.

Materials Provided:

Storage & Stability of Components:

1. FITC IgG H&L Conjugate No. AD CGEM2 (3.0 ml) with Evans Blue Counterstain is to be stored at 2-8°C upon receipt. The conjugate is stable at this temperature until expiration date on the vial label.
2. The antigen slides of M. Ovary sections must be stored at 2-8°C upon receipt. Check label for specific expiration date.
3. Universal negative control No. AD NC (1.0 ml) should be stored at 2-8°C upon receipt. Check label for specific expiration date.
4. Buffer Pack No. AD PBS1 - Phosphate Buffered Saline is stable at room temperature storage to the expiration date indicated on the label. The reconstituted Buffer does not contain preservatives and should be stored at 2-8°C. Care should be taken to avoid contamination.
5. Mounting Medium No. AD TMM3 is stable when stored at 2-8°C. Check label for specific expiration date.

Note: All kit components are available separately. Please see the current ALPHADIA Catalog for more details.

Additional Materials Required but not Provided:

Test tubes and rack or microtiter system
Disposable pipettes
Staining Dish and Slide Forceps
Moisture Chamber
Volumetric Flask (500 ml)
Distilled H₂O
Fluorescence Microscope
Paper Towels - lint free

Reagent Preparation:

1. Buffer Pack No. AD PBS1. Rehydrate buffer with 1 liter of sterile distilled water.

Specimen Collection:

Serological specimens should be collected under aseptic conditions. Hemolysis is avoided through prompt separation of the serum from the clot. Serum should be stored at 2-8°C if it is to be analyzed within a few days. Serum may be held for 3 to 6 months by storage at -20°C or lower. Lipemic and strongly hemolytic serum should be avoided. When specimens are shipped at ambient temperatures, addition of a preservative such as 0.095% sodium azide is strongly recommended.

Test Instruction:

Screening: dilute test serums 1:4 in PBS. **Titration:** set up doubling dilutions of serum starting at 1:4 1:8, 1:16, 1:32, etc.)

1. Once slides reach room temperature tear slide envelope at notch. Carefully remove the slide and avoid touching the antigen areas. The slide is now ready to use.
2. Place a drop of diluted serum (20 to 30 µl) and controls over the antigen wells.
3. Place slide with patient's serum and controls in a moist chamber for 30 minutes at room temperature (approximately 24°C).
4. Remove slide from moisture chamber and allow the serum to run off onto a piece of paper towel. Using a wash bottle, gently rinse remaining sera from slide being careful not to aim the rinse stream directly on to the well.
5. Wash in PBS for five minutes. Repeat using fresh PBS.
6. Place a blotter on the lab table with absorbent side up. Remove slide from PBS and invert so that tissue side faces absorbent side of blotter. Line up wells to blotter holes. Place slide on top of blotter. **Do not allow tissue to dry.** Wipe back of slide with dry lint free paper towel. Apply sufficient pressure to slide while wiping to absorb buffer.
7. Deliver 1 drop (25-30 µl) of conjugate per antigen well. Repeat steps 3-6.
8. Place 4-5 drops of mounting medium on slide.
9. Apply a 22 x 70 mm coverslip. Examine the slide under a fluorescent microscope.

Note: To maintain fluorescence, store mounted slide in a moisture chamber placed in a dark refrigerator.


Quality Control:

1. Positive and negative serum controls must be included in each day's testing to confirm reproducibility, sensitivity and specificity of the test procedure.
2. The negative serum control should result in little (1+) or no fluorescence. If this control shows bright fluorescence, either the control, antigen, conjugate or technique may be at fault.
3. The positive serum control should result in bright 3+ to 4+ fluorescence. If this control shows little or no fluorescence, either the control, antigen, conjugate or technique may be at fault.
4. In addition to positive and negative serum controls, a PBS control should be run to establish that the conjugate is free from nonspecific staining of the antigen substrate. If

the antigen shows bright fluorescence in the PBS control repeat using fresh conjugate. If the antigen still fluoresces, either the conjugate or antigen may be at fault.

Precautions:

1. All human components have been tested for (HbsAg) and HTLVIII/LAV by an FDA approved method and found to be negative. (Not repeatedly reactive). However, this does not assure the absence of HBsAg or HTLVIII/LAV. All human components should be handled with appropriate care.
2. The sodium azide (0.095%) included in the controls and conjugate is toxic if ingested.
3. Do not use components beyond their expiration date.
4. Follow the procedural instructions exactly as they appear in this insert to insure valid results.
5. For In Vitro Diagnostic Use.
6. Handle slides by the edges since direct pressure on the antigen wells may damage the antigen.
7. Once the procedure has started do not allow the antigen in the wells to dry out. This may result in false negative test results, or unnecessary artifacts.
8. Reagents should be inspected for evidence of bacterial or fungal contamination.
9. Use separate pipette tips for each sample and reagent to avoid cross contamination.
10. Do not reuse slides.

Component	AD PBS1 PBS Powder Packets AD TMM3 Mounting Medium	Precautionary Statement Prevention: P264 Wash thoroughly after handling. P280 Wear protective gloves and clothing.
Pictogram		Response: P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if easy to do. Continue rinsing.
Signal Word	WARNING	
Hazard Statement	H319 Causes serious eye irritation	P337+P313 If irritation persists, get medical advice/attention.

Bibliography:









1. Coulam, C.B. Prevalence of circulating antibodies directed toward ovaries among women with premature ovarian failure. Am J Reprod Immunol Microbiol 1985; 9: 23-4.
2. Mignot, M.H., Schoemaker, J., Kleingeld, M., Remanath, Rao B., Drexhage, H.A. Premature ovarian failure. I. the association with autoimmunity. Eur J Obstet Gynecol Reprod Biol 1989; 30: 59-66.
3. Gober, B., Barabarino-Monnier, P., Guillet-Rosso, F., Bene M.C., Faure, G.C. Ovary antibodies after IVF. Lancet 1990; 335: 723.
4. Coulam, C.B. The prevalence of autoimmune disorders among patients with primary ovarian failure. Am J Reprod Immunol 1983; 4:63-6.

ALPHADIA sa/nv
DIAGNOSTIC PRODUCTS
Avenue Vésale 26
B1300 WAVRE
BELGIUM
TEL : 32 (0) 10 24 26 49
FAX : 32 (0) 10 24 55 99
contact@alphadia.be



EC	REP	Emergo Europe Princessegracht 20 2514 AP The Hague The Netherlands
-----------	------------	---



	Manufactured by Prodotto da Fabricado por Fabriqué par hergestellt von
REF	Catalog number Numero di catalogo Número de Catálogo Numéro de catalogue Katalognummer
LOT	Lot Lotto Lote Lot Charge
EC REP	EC Authorized Representative Rappresentante Autorizzato CE Representante Autorizado CE CE Représentant autorisé EG autorisierter Bevollmächtigter
	EC Declaration of Conformity Dichiarazione di Conformità CE Declaración de Conformidad CE CE Déclaration de Conformité EG Konformitätserklärung
	Number of tests Numero di test Número de determinaciones Nombre de tests Anzahl der Tests
	See instructions for use Vedere le istruzioni per l'uso Consultar la instrucciones de uso Voir instructions Gebrauchsanweisung beachten
	Expiration date Data di scadenza Caducidad Date d'Expiration Halbbarkeitsdatum
 +2°C	Store at 2-8°C / 35-46°F Conservare a 2-8°C/35-46°F Almacenar a 2-8°C / 35-46°F Conserver à 2-8°C Bei 2-8°C / 35-46°F lagern
	Caution Attenzione Precaución Précautions Achtung
	Potential biological risk Potenziale rischio biologico Riesgo potencial biológico Biohazard Risque Biologique Potentiel Potentielle biologische Gefährdung
RFU	Ready for use Pronto all'uso Listo para su uso Prêt à l'emploi Gebrauchsfertig
IVD	For in vitro diagnostic use Per uso diagnostico <i>in vitro</i> Para uso solo in vitro Usage in vitro Für in-vitro diagnostische Verwendung
RUO	For research use only Solo per ricerca Para uso solo en investigación Pour recherche Nur für Forschungszwecke
IUO	For investigational use only Solo per uso investigativo Para uso solo en investigación Pour investigation Nur für Forschungszwecke
IFA/DFA PBS	Phosphate Buffered Saline Tampone salino fosfato Fosfato Salino Tamponado Tampon phosphate salin PBS
SOR	Sorbent Assorbent Sorbente Absorbant Sorbens

SLIDE	Tissue Substrate Slide Vetrini con substrate di tessuto Porto objetos de Substrato de Tejido Lame portant le substrat tissulaire Gewebe substrat-Objektträger
MM	Mounting Medium Mezzo di montaggio Medio de Montaje Liquide de montage Eindeckmedium
10x	Concentration Concentrazione Concentración Concentration Konzentration
ENS	Enhancement solution Soluzione di rinforzo Solución de realce Solution d'amplification Verstärkungslösung
WASHB	Wash Buffer Tampone di lavaggio Tampón de lavado Tampon de lavage Waschpuffer
MPS 12x8	Microplate Strips Strip per Micropiastra Tiras de micro placa Microplaque Mikrotiterplattenstreifen
CONJ	Conjugate Coniugato Conjugado Conjugué Konjugat
SUB	Substrate Substrato Sustrato Substrat Substrat
STOP	Stop Solution Soluzione bloccante Solución de Parada Solution d'arrêt Stopplösung
CAL X	Calibrator(s) Calibratore (i) Calibrador (s) Calibrateur(s) Kalibrator(en)
CONTROL -	Negative Control Controllo Negativo Control Negativo Contrôle Négatif Negative Kontrolle
CONTROL +	Positive Control Controllo Positivo Control Positivo Contrôle Positif Positive Kontrolle
CONJ CNS	Counterstain Colorante di contrasto Contraste Contre colorant Gegenfärbung
CS	Coverslip Coprioggetto Cubre portaobjetos Lamelles couvre-objet Deckgläser
CONJ +	Positive Conjugate Coniugato Positivo Conjugado Positivo Conjugué Positif Positivekonjugat
CONJ -	Negative Conjugate Coniugato Negativo Conjugado Negativo Conjugué Négatif Negativkonjugat
DIL	Sample Diluent Diluente del campione Diluyente de muestra Tampon de dilution Probenverdünnungslösung